



True Color Texture, Rpd Expert

Highly Restore The Real Color Of The Model To Help

The Rpd Design

Texture Expert To Simplify The Rpd Process

Provide One-stop Solution For Digitalization Of Rpd



Model	UP560
Size	305*310*550mm
Weight	10kg
Cameras	2*3.0mp
Accuracy	6 microns
Projector	Customized blue-light / Multi-color scanning
Interface	USB 3.0
Output Formats	STL, UM, PLY, OBJ
Power Supply	AC 110-220V, 50Hz
Product Features	All-in-one scan, One click scanning triple tray
	Detect margin line, Texture scan
Scan Speed	Full arch: 12s, Unsectioned models scan: 113s
	post&core scan: 93s, All-in-one: 11s
	Triple tray: 101s, True color texture scan: 126s



