## **UPCERA**



**C41** 

4-Axis Wet Milling Machine

Reliable Materials Expert



## C41

## 4-Axis Wet milling Machine

- Compact and ergonomic design
- Built-in automatic circulating water system
- Rapid and accurate production with added efficiency
- Versatile indications up to 40 mm block
- Ideal for lab and clinic for same-day restorations



## **Device Paremeters**

Specification	Parameters
Dimensions	45*35*42cm
Weight	45kg
Voltage	220V 50Hz
Power	<1KW
Compatible Compressed Air	0.42-0.8 MPa
Spindle Speed	80,000rpm
Precision	< 0.01mm
Number of Bur	5
Cooling Type	Water cooling
Materials	Composite Resin Block, Glass ceramic, Titanium Block PMMA, Hyramic
Application	Veneer, full contour crown, inlay, onlay, pre-milled abutment
Milling Speed	PMMA (12 minutes/pc), Glass-ceramic (10~30 minutes/pc), Premill (40 minutes/pc)