UPCERA



FLNT Pearl

Zirconia Sintering Beads

Reliable Materials Expert



FLNT Pearl

Zirconia Sintering Beads

FLNTPearl can withstand up to 1700°C and are made of high-purity zirconia that provide contamination-free results for restorations.



Homogeneous size and surface quality

FLNTPearl Zirconia Sintering Beads provide a better temperature rate when using sintering trays. Because the sintering tray and FLNTPearl heat at two different temperature rates, FLNTPearl absorbs the heat and provides a controlled temperature which helps prevent deformation. Zirconia Sintering Beads support crowns, inlays, onlays, veneers, and bridges during the sintering process.