

## **Grinding and Polishing kits**

Reliable Materials Expert

## Ceramic Diamond Grinders Special for Zirconia

Material: Diamond grits mixing with ceramic powder.

Usage: For zirconia inner crown fast finishing, dry grinding, without water cooling (it is a effective fast way to instead of the diamond burs with water cooling way)

Caution: Please use the products in allowed RPM.OR it may take the products damaged!

Grit: Mid-Coarse Packing(pcs): 1	Ţ		-	9	Ţ	Ť	+
Order No.:	CD 001 3	CD 006 3	CD 008 3	CD 007 3	CD 009 3	CD 011 3	CD 012 3
Recommended R.P.M:	15,000-30,000 / min	15,000-30,000 / min	15,000-30,000 / min	15,000-30,000 / min	15,000-30,000 / min	15,000-30,000 / min	15,000-30,000 / min
Size Φ × L(mm):	13 x 7	16 x 2	22 x 2	16 x 3.5	22 x 4.5	16 x 1.5	22 x 2
ISO No.: 805	900 030 524 130	900 372 524 160	900 372 524 220	900 041 524 160	900 041 524 220	900 307 524 130	900 307 524 220
Grit: Mid-Coarse Packing(pcs): 1						Ĭ	
Order No.:	CD 201 3	CD 202 3	CD 203 3	CD 212 3	CD 214 3	CD 221 3	CD 222 3
Recommended R.P.M:	20,000-35,000 / min	20,000-35,000 / min	20,000-35,000 / min	20,000-35,000 / min	20,000-35,000 / min	20,000-35,000 / min	20,000-35,000 / min
<u>Size Φ × L(mm):</u> ISO No.: 805	3 x 7 104 110 524 030	5 x 12 104 113 524 050	7 x 12 104 113 524 070	4 x 10.5 104 174 524 040	3 x 7 104 197 524 030	5 x 2.5 104 014 524 050	13 x 3.5 104 024 524 130
Grit: Mid-Coarse Packing(pcs): 1		Ĩ	-	-	ł		
Order No.:	CD 224 3	CD 225 3	CD 231 3	CD 232 3	CD 251 3		
Recommended R.P.M:	20,000-35,000 / min 3.6 x 7	20,000-35,000 / min 6 x 8	20,000-35,000 / min 12 x 2	20,000-35,000 / min 12 x 3	20,000-35,000 / min 3 x 3		
Size Φ × L(mm): ISO No.: 805	104 237 524 036	104 030 524 060	104 041 524 120	104 042 524 120	104 001 524 030		
Packing(pcs): 1 Order No.: Recommended R.P.M: Size Φ × L(mm): ISO No.: 805	CD 232 4 20,000-35,000 / min 12 x 3 104 042 514 120	CD 202 4 20,000-35,000 / min 5 x 12 104 113 514 050	CD 212 4 20,000-35,000 / min 4 x 10.5 104 174 514 040	CD 214 4 20,000-35,000 / min 3 x 7 104 174 514 030			

## Rubber Diamond Polishers Special for Zirconia/Porcelain Workpiece

Material: Diamond grits mixing with synthetic rubber

Usage: For smoothing and polishing zirconia inner crown, zirconia veneer/porcelain teeth fast super bright. (Avoid the color differences and porcelain fragile caused by second glaze and baking. And it also can save time)

Caution: Please use the products in allowed RPM.OR it may take the products damaged!

Grit: Packing(pcs): 1	Extra Coarse	Mid-Coarse	Fine	Extra Fine	Ultra Fine
Order No.:	RD 011 1	RD 011 3	RD 011 4	RD 011 5	RD 011 6
Recommended R.P.M:	15,000-25,000 / min	15,000-25,000 / min	15,000-25,000 / min	10,000-20,000 / min	10,000-20,000 / min
Size Φ × L(mm):	19 x 2	19 x 2	19 x 2	19 x 2	19 x 2
ISO No.: 802	900 307 544 190	900 307 524 190	900 307 514 190	900 307 504 190	900 307 494 190
Grit: Packing(pcs): 1 Order No.:	Extra Coarse RD 012 1	Mid-Coarse RD 012 3	Fine RD 012 4	Extra Fine RD 012 5	Ultra Fine RD 012 6
Recommended R.P.M:	15,000-25,000 / min	15,000-25,000 / min	15,000-25,000 / min	10,000-20,000 / min	10,000-20,000 / min
Size $\Phi \times L(mm)$ :	26 x 3	26 x 3	26 x 3	26 x 3	26 x 3
ISO No.: 802	900 307 544 260	900 307 524 260	900 307 514 260	900 307 504 260	900 307 494 260

Grit: Mid-Coarse Packing(pcs): 1	4	F		t	t	t	t
Order No.:	RD 201 3	RD 202 3	RD 211 3	RD 212 3	RD 216 3	RD 221 3	RD 222 3
Recommended R.P.M:	15,000-25,000 / min						
Size Φ × L(mm):	3 x 7	5 x 12	3 x 7	4 x 12	4 x 12	9 x 8.5	18 x 3
ISO No.: 802	104 110 524 030	104 113 524 050	104 243 524 030	104 243 524 040	104 174 524 040	104 030 524 090	104 026 524 180

Grit: Mid-Coarse Packing(pcs): 1	T	T	Ţ	Ŧ	Ţ
Order No.:	RD 231 3	RD 232 3	RD 233 3	RD 241 3	RD 237 3A
Recommended R.P.M:	15,000-25,000 / min				
Size Φ × L(mm):	12 x 1	12 x 2	12 x 4	12 x 2	15×2.5
ISO No.: 802	104 371 524 120	104 041 524 120	104 042 524 120	104 304 524 120	104 041 524 150